










Edgeコーティング電極 - ディスポーザブル (絶縁付) -

Edge Coated Electrodes

製品コード		全長	通電部分	発注単位	クラス
E1455		7.0cm	5.1mm	50本/箱	II
E1455B		7.0cm	5.1mm	25本/箱	II
E1455-4		10.16cm	5.1mm	50本/箱	II
E1455B-4		10.16cm	5.1mm	25本/箱	II
E1455-6		16.51cm	5.1mm	50本/箱	II
E1465		7.2cm	5.1mm	50本/箱	II
E1465B		7.2cm	5.1mm	25本/箱	II
E1465-4		10.16cm	5.1mm	50本/箱	II
E1465B-4		10.16cm	5.1mm	25本/箱	II
E1465-6		16.51cm	5.1mm	50本/箱	II

実寸大

Edgeコーティング電極 – ディスポーザブル –

Edge Coated Electrodes

製品コード		全長	通電部分	発注単位	クラス
E1450X		6.35cm	2.8cm	50本/箱	II
E1450G		7.62cm	2.8cm	50本/箱	II
E1475X		7.0cm	2.8cm	50本/箱	II
E1452		7.2cm	2.8cm	50本/箱	II
E1450-4		10.16cm	2.54cm	50本/箱	II
E1450-6		16.51cm	2.54cm	50本/箱	II
E1452-6		16.51cm	2.54cm	50本/箱	II

実寸大

Medtronic

コヴィディエンジャパン株式会社

Tel: 0120-998-971

medtronic.co.jp

販売名：フォース電気手術器
医療機器認証番号：221AABZX00018000

ct-he-ea-edge(01)
1711.pdf.cry